Covid-19 Testing CareEvolve Portal

CRSP

CRSP CareEvolve Website


Employer/School - Facing Website

- Manage Participant roster
- Manage frequency of testing
- See Results

Participant - Facing Website


Access your results

Employer Facing Website



Employers - Load Participant Roster

1 - .CSV upload




2 - Manual Entry



Required Fields:

- First Name
- Last Name
- Date of Birth
- Gender
- Address
- Phone Number
- Email Address
- Race
 - Ethnic group


Employers - Create Orders (Set frequency for testing)




Define the Start and End dates of a Standing Order, and then choose the frequency.

After the schedule summary page appears, the individual orders will be listed to view or change.

Employers - Create Orders (Set frequency for testing)




Employers - Search for Employees, Complete Orders when ready to collect




Employers - Search for Employees, Complete Orders when ready to collect

Complete an Order


Click the blue [Draw Now] button.



Click the [Complete] button.




Employers - Print Label or Order Requisition





Collection 'kit' (provided to your site)




Ship rack of tubes to **Broad CRSP**

Testing


Collect sample

Employers - See Results within 24 hours









SARS-CoV2 Real-time Reverse Transcriptase (RT)-PCR Diagnostic Assay Positive Negative C -- END OF REPORT --

© CareEvolveTM Result Activities

< Back to Results Inbox

Chart Tools **Print Result**

Fax -

Email -

Result 1 of 5

Employers - Search all results and view all positive / negative cases

Global Results Search Back to Global Results Search Lab: CRSP Location: All Locations Age: All Any analyte(s) Below: Covid19_Diagnostic [SARS-CoV2 Real-time Reverse Transcriptase (RT)-PCR Diagnostic Assay]: Abnormal

11002

Name	Requisition	Office Id	Birth Date	Phone	Provider	Test	Result	Coll Date	
Test, Kim	CQX6U6BPFP11		01/25/1989	988-888- 8888	Test, Physician	SARS-CoV2 Real-time Reverse Transcriptase (RT)-PCR Diagnostic Assay	Positive	05/23/2020 16:06	0
Test, Test	CQX6U6BPFP1	OFFICEID1234	01/01/1980	888-888- 8888	Test, Physician	SARS-CoV2 Real-time Reverse Transcriptase (RT)-PCR Diagnostic Assay	Positive	05/23/2020 16:06	0
Test, Test	CQX6U6BPFP1RS0G	OFFICEID1234	01/01/1980	888-888- 8888	Test, Physician	SARS-CoV2 Real-time Reverse Transcriptase (RT)-PCR Diagnostic Assay	Positive	05/23/2020 16:06	0
Test1, Kim1	CQX6U6BPFP111		01/26/1989		Test, Physician	SARS-CoV2 Real-time Reverse Transcriptase (RT)-PCR Diagnostic Assay	Positive	05/23/2020 16:06	0
1	•								

Participant Facing Website

Participants will be emailed when their results are ready

Your recent laboratory results are available from your provider, please register now.

Registration Code: CE717D5302

Please use the code above and the following link to register:

https://crsp.careevolve.com/Patient/Profile/RegStart

If you have any questions regarding this email, please contact Customer Support.

https://crsp.careevolve.com/Patient/Home/Contact

